



MicrodermMD
Home
Microdermabrasion
System



CE



Congratulations, you have taken the first step towards a dermatologist-grade skincare regimen that is fast, easy, and effective.

Triple Action System

Use the powerful Vacuum Suction Motor, natural Diamond Tip, and Blackhead Extraction Tip to get professional results at home.

Get ready for the skin you've always dreamed of, using the safe and effective MicrodermMD, you can now: Improve skin tone and texture smooth the appearance of fine lines. Reduce the appearance of age spots or sun damage. Have younger, more radiant looking skin. This dermatologist recommended treatment is the skin maintenance tool you need to keep your skin looking young and refreshed.

Read on for the simple instructions on how to use your new MicrodermMD device.

MICRODERM MD DEVICE OVERVIEW



PRECAUTIONS

1. Do not touch the power outlet or adaptor plug with wet hands
2. Do not disassemble the MicrodermMD for any reason, contact us if you need assistance
3. Keep the custom MicrodermMD black wool filter in at all times to prevent foreign substances from entering the suction motor
4. Do not drop the device or cause any high impact as it may damage the frame or motor
5. Be careful not to excessively press the Diamond Tip to your skin during usage, start with one or two passes and increase usage over time
6. Do not use the diamond tip around the eye area, any areas that contain inflammatory acne, or any areas that have skin rashes or wounds
7. Soothe your skin by washing with cold water and apply a light moisturizer after each treatment
8. Stop using the device immediately if you develop any abnormal reactions such as itching or bothersome inflammation
9. Make sure to apply sunblock after treatments if you will be exposed to the sun
10. Do not use the MicrodermMD device for commercial purposes, it is intended for personal use at home

SPECIFICATIONS

- What's included: MicrodermMD unit, Diamond Tip, Blackhead Extraction Tip, cleaning brush, 100 pack of custom MicrodermMD wool filters, instruction manual, power adaptor
- Specification: 300 X 230 x 75
- Weight: 1.3 kg/unit
- Input power: AC Single Phase 100-220V
- Frequency: 60 Hz
- Power consumption: 13 Watts



PRE-EXISTING CONDITIONS AND PRECAUTIONS

Do NOT use the MicrodermMD without first consulting a skincare professional if you have any of the following pre-existing conditions:

Lesions, warts, or recent herpes outbreak

Sunburn or sun sores

Active rosacea

Auto-immune disorders

Skin cancer

Vascular lesions

Any open wounds, sores, or cuts

Taking oral blood thinners (anti-coagulants)

HOW TO USE

Step 1.

Plug the power adaptor into the back of the MicrodermMD unit, then plug the adaptor into the wall. The red light above the "Power" button should now light up.

Step 2.

Press the "Start/Stop" button to start the suction motor. The default setting is Suction Level 1. We recommend starting with Suction Level 1 or 2 for the first few uses to allow your skin to get accustomed to the microdermabrasion treatment.

Step 3.

With the suction motor running, place the diamond tip against the area to be treated. It can help to pull the skin tight with your free hand, especially in hard to treat areas like the nose. The diamond tip must be FLAT against the skin for proper suction to take place.

Step 4.

Stroke the diamond tip against the skin with smooth, even motions, making sure not to stroke more than 3 times in any given area.

Step 5.

Treat areas such as the face, neck, chest, arms, hands, feet, or any other body area 2-3 times per week using this method. You can increase the Suction Level over time to Suction Level 5 as long as you are within your comfort level.

Step 6.

On days that you have NOT used the Diamond Tip to treat an area, you may instead use the Blackhead Extraction Tip. This tip creates a very powerful suction that extracts the oil, sebum, and buildup inside the pores. It can be used on any area of the body, and is especially useful on the top and sides of the nose.

EXAMPLE TREATMENT PROTOCOL FOR FACE

1. Cheeks (Level III)
 - a. Move the diamond tip a short distance outward from the nose towards the sides of the cheeks
 - b. Move the diamond tip in a downward motion towards the lips
2. Around the mouth (Level II)
 - a. Move the diamond tip downward over the upper lip and outward around the lips
 - b. Cover the chin area in downward strokes and around the lower lip region
3. Around the eyes (Level I)
 - a. Lightly move the tip around the eye region with outward strokes only
 - b. Move the tip along the ends of both eyebrows
4. Forehead (Level IV)
 - a. Move up the center of the forehead and outwards towards the temples
 - b. Use up and down motions along the width of the forehead
5. Around the nose (Level IV)
 - a. Move down along the ridge of the nose, making sure to maintain suction at all times
 - b. Move the tip sideways across each side of the nose

SKIN CARE AFTER TREATMENT

You may experience some redness and sensitive skin after your treatments, this is completely normal. This will usually subside within a few hours, and will vary greatly from person to person. You can lower the Suction Level or the amount of weekly treatments to suit your own comfort level. Microdermabrasion is a gradual treatment process that will continue to improve the appearance and health over your skin for years.

Step 1. Wash your face with cold water to remove dead skin cells and debris.

Step 2. Apply moisturizer or sunblock after each treatment.

Step 3. Make sure to remove the black wool filter and replace it with a new one after each treatment.

The Diamond Tip that comes with your MicrodermMD unit is made of real diamond powder, making it very durable. The Diamond Tip is built to last at least 1 year with normal usage, however, it will need to be cleaned with an Ultrasonic Cleaner every few months. You can purchase an Ultrasonic Cleaner very inexpensively online, please contact us for details. Also, many jewelry stores will do the cleaning at no charge or for a small fee. We highly recommend this cleaning as it will keep the Diamond Tip in new condition.

CARING FOR THE DEVICE

- Keep the device away from water and humidity
- Be careful to prevent foreign substances from entering the motor by keeping a black wool filter in at all times
- Clean the diamond tip with an Ultrasonic Cleaner every few months
- Replace the black wool filter after every treatment
- Be sure to unplug the power cord when the MicrodermMD is not in use

TROUBLESHOOTING TIPS

Issue: MicrodermMD does not turn on

Solution: Make sure the power adaptor is plugged in securely at the back of the unit, and plugged into the wall power outlet. Check to see if the light is on above the "Power" button. If the light does not come on, please contact us for a replacement power adaptor.

Issue: MicrodermMD has lost suction power

Solution: Check the Diamond Tip closely for any hairline cracks, punctures, or spaces that would cause a loss of suction. Remove the Diamond Tip and ensure that both rubber O-rings are in place. If either problem exists, contact us for a replacement Diamond Tip or new O-rings.

Issue: Diamond Tip is no longer coarse enough

Solution: The Diamond Tip will accumulate skin cells in between the diamond grains over time, making the tip feel smoother and less effective at exfoliating the skin. We recommend using an Ultrasonic Cleaner to remove the skin cells every few months. You can purchase an Ultrasonic Cleaner for use at home, or take the diamond tip to a jewelry store for professional cleaning. This will take the Diamond Tip back to new condition.

SATISFACTION GUARANTEE

Trophy Skin is committed to 100% Customer Satisfaction, and unlike most companies, we put our money where our mouth is. If for any reason you are unsatisfied with the results of your MicrodermMD device, you may initiate the return process by sending an email to service@trophyskin.com within 30 days of receiving your device. The following terms will apply to be eligible for our full, "no questions asked" refund policy:

- MicrodermMD device is in new condition
- Device must be shipped back, at your cost, using the shipping method of your choice, as long as a tracking number is provided
- MicrodermMD must be returned in the original box, with all included accessories, power supply, and original packing materials



Clearer skin provided through microdermatation

WARRANTY

Your MicrodermMD device is warranted to be free from defect in material and workmanship for a period of three years under normal use after its original purchase date. This warranty extends only to the original retail purchaser and only when purchased through the Trophy Skin website or approved 3rd party affiliated sites.

If your MicrodermMD device should become defective during the warranty period, contact Trophy Skin via email at service@trophyskin.com for repair or replacement. We reserve the right to replace a defective product with the most comparable product currently available. We will work with you to ensure that you are taken care of to the best of our ability, our first priority is to make sure you have a working device that you can use.

If you have any warranty issues, claims, or questions, please give us a call at 469-893-0710 and we will resolve the issue for you promptly.

DISCLAIMER

Trophy Skin Is Not Engaged in the Practice of Medicine

Trophy Skin expressly disclaims any and all responsibility for any consequences, directly or indirectly, related to any use, misuse or interpretation of the web site or any other information provided or omitted in connection therewith and you agree to indemnify, defend and hold Trophy Skin harmless therefrom. Trophy Skin shall not be deemed to be engaged in the corporate practice of medicine as a consequence of selling the MicrodermMD. In addition, I specifically agree to with the following statements:

- I understand the MicrodermMD is a cosmetic treatment and that no medical claims are expressed or implied.
- I understand that to achieve maximum and continued results the protocol recommended by Trophy Skin should be followed.
- I understand there are no guarantees implied as to the results of MicrodermMD treatment, due to many variables, such as: age, skin type, skin condition, sun damage, smoking, alcohol, environmental exposures, etc. I understand that I may or may not actually see demonstrable visible results, that each case is individual.
- I understand that I will not be allowed to have treatments during any pregnancy.

MEDICAL TREATMENT OR DIAGNOSIS DISCLAIMER:

Nothing contained in www.trophyskin.com is intended to be instructional for medical diagnosis or treatment. The information presented by www.trophyskin.com and/or Trophy Skin should NOT be considered complete, nor should it or take the place of a consultation with a physician or competent healthcare professional for medical diagnosis and/or treatment, but is provided as information only.

Information obtained on this site related to anti-aging treatment is NOT meant to be exhaustive and obviously can not cover all situations and conditions. Again, we urge you to seek consultation or advice of your physician or other qualified health care provider promptly with any health care related questions. You should never disregard medical advice or delay in seeking it because of something you have read on this web site or others.

The information contained in www.trophyskin.com is compiled from a variety of information providers. The web site www.trophyskin.com does not directly or indirectly practice medicine or dispense medical services.



 **Trophy**
Skin